



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Binocular Stereo Microscope

**Product Code :**  
SAMCOMEDICAL4011A003



**Description :**

Samcomedical Binocular Stereo Microscope

**Technical Specification :**

Binocular Stereo Microscope Binocular Stereo Microscope a significant advantage of Stereoscopic Microscopic is that by using two optical systems, an erect image of the object can be seen in three dimensions. Magnifications upto 80x are possible by using additional eyepieces. Stereo Microscope are used in every field of Education, Electronics, Watch Industry, Tool Industry, Food Industry, Dissection of Organisms, Medical Science, Quality Control, Inspection & Research Purposes. Observation Head - Binocular head inclined at 45-deg., rotated through 360-deg with adjustable inter-pupillary distance. Focussing - By Rack and Pinion mechanism with adjustable tension collar. Magnification - Standard magnification 20X & 40X or 10x & 30x. Illumination - With integral illuminator, 2 lamps are fitted to provide both transmitted and incident illumination. Objectives - 2 sets of objectives 2X and 4X, selected by rotating nose piece through 90-degrees. Eyepieces - Standard extra wide field eyepieces WF-10X is provided. Binocular Stereo Microscope AEI-C, Binocular Loupes exporters, Binocular Loupes manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

