



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fujitron Zeiss Type Slit Lamp

Product Code :
SAMCOMEDICAL4002A042



Description :

Fujitron Zeiss Type Slit Lamp Fujitron Zeiss type Slit lamp with high-grade multi-coated optics and its elegant design realizes high transmittance with excellent resolution image without ghost and flare The classical Galilean Converging Binocular tube with high eye point eyepieces provides wide field of view The Drum style 3step & 5step magnification control can achieve from low to high magnification as required by the Surgeon, eliminating the need of additional Eyepiece Various kinds of Built-in filters like Cobalt blue (excitation), Green(Red-free), Infrared cut(Heat absorbing), filters can be selected through a filter lever High luminance Halogen lamp, with intensity control knob on the slit lamp base (optional) enhances the easy accessibility while using All Fujitron slit lamp are designed to meet the required specifications for future up gradation to Digital imaging or video Slit lamps. Fujitron Zeiss Type Slit Lamp, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,