



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Automatic Capsule Filling Machine

Product Code :
SAMCOMEDICAL4003A001

Description :

Automatic Capsule Filling Machine We offer a precision engineered range of automatic capsule filling machines that are skillfully designed to meet the precision production demands of modern pharmaceutical procedures. These machines are ideal for filling hard gelatin capsules of various sizes with powder, granules and pellets. This machine is productivity enhancing machine which is handled by one operated, our range is also available in customized specifications to the clients. Salients features: Compact design to cGMP Standard Filled weight accuracy + 2% Total PLC Controlled Machine Varied output range of 6,000,25,000,40,000,90,000,150,000 capsules per hour.Total PLC Controlled Machine Quick & easy change over of capsule size from '0' to '4' Individual air control for faulty & filled capsulesAutomatic Capsule Filling Machine, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,