



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
TABLET FRIABILITY TEST MACHINE

**Product Code :**  
SAMCOMEDICAL4003A018



**Description :**

TABLET FRIABILITY TEST MACHINE This machine is suitable to find out the rate of abrasion and impact hardness of tablets. The unit is electrically operated and is supplied complete with timer, 12" dia. Plexiglass drum etc TABLET FRIABILITY TEST MACHINE, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),