



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
P.P. CAPPING / T.O. OR T.D. CAPPING MACHINE [SEMI-AUTOMATIC

Product Code :
SAMCOMEDICAL4003A027



Description :

P.P. CAPPING / T.O. OR T.D. CAPPING MACHINE [SEMI-AUTOMATIC American type industrial model (motorized light duty) incorporating rotary threading with sealing rollers and collect head complete with 3 sealing heads for 3 different- sizes of caps, ½ HP Motor, M.S. Stand claddcd with S.S. with S.S. Top, Acrylic Shield, Belt Guard, (Made out of S.S.) On-Off switch & Paddle Lifting arrangement. The machine is as per -GMP Standard, No painted parts visible.



Samco Medical

Website: www.samcomedical.com,