



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
DIAMOND WHEELS

Product Code :
SAMCOMEDICAL4004A040



Description :

DIAMOND WHEELS For long life reverse / face direction of the wheel after 1 to 2 month The hole / bore of the wheel should fit in bearing position of the shaft Regular use of proper dressing stick on wheel increases cutting rate and also increases life of wheel DIAMOND WHEELS, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,