



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Delivery Bed With Removable Legs Section

**Product Code :**  
SAMCOMEDICAL10014A003

**Description :**

Delivery Bed With Removable Legs Section

**Technical Specification :**

Delivery Bed With Removable Legs Section Delivery Table / Bed, made in 3 sections. Body Section with legs mounted on rubber tips, legs section on swivel castors with brakes. Size: 1680 i<sub>l</sub> 1/2 750 i<sub>l</sub> 1/2 800 mm. With adjustable back rest to varying positions. Framework made of Mild Steel Epoxy coated finish. Supplied with S.S. bowl and two leg holders. Mounted on protective stumps. Freight saving knock-down construction.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),