



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Emergency Two Side Trolley

**Product Code :**  
SAMCOMEDICAL10015A018

**Description :**

Emergency Two Side Trolley

**Technical Specification :**

Emergency Trolley Two side Trolley designed to supply medicines. (Front & back same) Size: 850  $\frac{1}{2}$  515  $\frac{1}{2}$  985 mm. 4 drawers with dividers suitably designed for keeping medicines. 1 Multi purpose drawers. Equipped with IV pole, waste bin, needle disposable container, file cassette & guard rails. Equipped with lock key system. 4 swivel castors (2 with brakes).



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),