



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Foot Stool (Single step)

Product Code :
SAMCOMEDICAL10024A004

Description :

Foot Stool (Single step)

Technical Specification :

Foot Stool (Single step) Size : 51L x 30W x 23H cms. ERW steel tubular frame work mounted on PVC shoes. Top fixed with antistatic rubber mat surrounded by aluminium angle frame. Pre-treated and epoxy powder coated.



Samco Medical

Website: www.samcomedical.com,