



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Nasal Oxygen Cannula

Product Code :
SAMCOMEDICAL500001A052



Description :

Nasal Oxygen Cannula Used for efficient administration of oxygen for greater comfort of patient. Twin Nasal Tips provide maximum freedom to the patient and leaves the mouth free. Nasal Tips ensures equal volume of oxygen to both air passage. Manufactured from soft non toxic, medical grade PVC Compound. Multichannel tube ensure regular supply of oxygen if tube is accidentally kinked. Proximal end of tube provided with connector for easy connection with oxygen source. Non-Sterile, individually packed in a polybag



Samco Medical

Website: www.samcomedical.com,