



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi-Fowler Bed

Product Code :
SAMCOMEDICAL1001A034

Description :

Semi-Fowler Bed

Technical Specification :

Semi-Fowler Bed Back rest section maneuvered smoothly by screw mechanism (hand crank) from foot end. Overall size : 2030L x 900W x 600H mm. Frame work of mild steel rectangular CRC tubes. Mounted on 10 cms. dia heavy duty castors (2 with brakes). Rigid mesh mattress platform. Stainless steel head and foot bow with laminated wood panels. Provision for saline / IV rod at four locations. Pre-treated and epoxy powder coated finish. Knockdown construction



Samco Medical

Website: www.samcomedical.com,