



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Overbed table, Height Adjustable

**Product Code :**  
SAMCOMEDICAL10014A033



**Description :**

Overbed table, Height Adjustable Overbed table with Laminated Wooden Top. Size: 795 ½ x 400 ½ x 710 / 970 mm. Framework is made of Epoxy steel. Adjustable height from 710 - 970 mm by gas spring. 4 Swivel castor



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),