



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
LED Shadowless OT Light Single Reflector Stand Floor
Model With Castors

Product Code :
SAMCOMEDICAL200001A077



Description :

LED Shadowless OT Light Single Reflector Stand Floor Model With Castors Illumination Intensity : 40000-65000 Lux Color Temperature : 3800-5000 K Color rendering index : 85 (Ra) Spot diameter : 80-160mm Height Adjustment : 700mm Temperature Rise of Field : 0 ($\frac{1}{2}^{\circ}\text{C}$) Specification of Bulb : LED Bulb 5x3W Power Input : 15 VA Power Supply : 220 V / 50 Hz. AC Brightness Control : Available



Samco Medical

Website: www.samcomedical.com,