



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Suction Unit Eurovac A Battery Operated

**Product Code :**  
SAMCOMEDICAL150001A002



**Description :**

Suction Unit Eurovac A Battery Operated

**Technical Specification :**

Suction Unit Eurovac A (Battery Operated) Housing: MS Powder Coated Capacity: -600 mm Hg 10 , 14 Ltrs/ min  
Pump Type: Piston Pump Jars: 1 x 0.5 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable  
Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60  
Hz, 55Watt , 12 V DC (Internal Rechargeable Battery 12V 4.5Amp) Battery: Lead Acid (Ni-cd Optional) Noise  
Level: 55 dB A 3 Dimension & Weight: 37 x 14 x 30 cms, 5.5 Kg. Overflow Safety: Mechanical / Electronic  
Acoustic Alarm Accessories: Car Charger & Power Adaptor



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),