



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
C-Arm Compatible Fully Electromatic O.T. Table With  
Accentric Base

**Product Code :**  
SAMCOMEDICAL9007A005



**Description :**

C-Arm Compatible Fully Electromatic O.T. Table With Accentric Base Technical Constraints Top Dimension : L 1930 xW 553mm Height Adjustmen : 762mm- 1012mm Trendelenburg / Reverse : 30° / 25° Lateral Tilt : 20° / 20° Kidney Elevator : 150mm Back Rest (up / down) : 80° / 25° Leg Rest (up / down) : 15° / 90° Head Rest (up / down) : 20° / 60° Power Supply : 24V DC Battery Backup (Optional at extra price) : 2 - 3 hours Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs Overview : Electromatic O.T. Table is equipped with imported electric actuators for smooth & effortless working of the table. It ensures the actuators working in a constant speed without any noise & is designed to run even at low voltage to carry on the surgery Key Structure: The table has been made from acid proof, stainless steel (the face of the table is easy to clean and immune to disinfection agents) Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses. The various positions such as Up - Down movements, Trendelenburg - Reverse Trendelenburg, Lateral Tilts, Flex - Reflex, Chair Position as realized by remote control



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),