



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Side End Contolled Universal Operating Table

**Product Code :**  
SAMCOMEDICAL9007A019



**Description :**

Side End Contolled Universal Operating Table Standard Features: Suitable for thoracic surgery, abdominal surgery, ophthalmology, otorhinolaryngology, obstetrics and gynaecology, urinary surgery and orthopaedics. Height adjustment is hydraulically operated and the rest of the movements are controlled mechanically. Accessories and framework of the table made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. The Leg section is dismountable and separable Technical specifications Table Top (Length x Width) : 2100mm x 480mm Height Adjustable : 800mm-1045mm Trendelenburg :  $\pm 35^{\circ}$  Reverse Trendelenbug :  $\pm 20^{\circ}$  Lateral Tilt :  $\pm 22^{\circ}$  /  $22^{\circ}$  Back Rest (up/down) :  $\pm 75^{\circ}$  /  $22^{\circ}$  Head Rest : It can be folded upto  $90^{\circ}$  Kidney Bridge Elevation : 120mm Leg Rest adjustment (down) :  $90^{\circ}$  Standard Accessories: Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),