



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Frog Splint

Product Code :
SAMCOMEDICAL100004A003

Description :

Frog Splint

Technical Specification :

Frog Splint Size: S, M, L Feature: Foam padded aluminium splint in unique shape to cover the dorsal surface of finger to remain open Four extended arms to ensure good fitting. Application: Hyper flexion injuries Swan neck deformity To hold distal inter phalangeal joints in correct alignment Collateral ligament injuries of distal inter phalangeal joint



Samco Medical

Website: www.samcomedical.com,