



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Surgical Endoscopic Ceiling Pendant

Product Code :
SAMCOMEDICAL4007A008



Description :

Surgical Endoscopic Ceiling Pendant Key features : Cross Arm and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Horizontal rotation can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate position locking. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal. High bearing capacity crane for suspending endoscopic dolly and cavity mirror Basic configuration : One Gas discharge terminal for compressed air, oxygen, carbon dioxide, suction and anesthetic gas Four power outlets



Samco Medical

Website: www.samcomedical.com,