



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Surgical Operation Light - Superior (4
Bulbs, 500mm Dome Dia.)

Product Code :
SAMCOMEDICAL4007A028



Description :

Ceiling Shadowless Surgical Operation Light - Superior (4 Bulbs, 500mm Dome Dia.) Light head : 500 mm dome dia., 4 bulbs. Light Intensity : 75,000 Lux. Halogen Bulb : 24V x 50W x 4 Nos. Power Consumption : 200VA. Action Radius : 1,200 mm. Colour Temperature : 4200 \pm 1/2 K. Temperature Rise : 5 \pm 1/2-8 \pm 1/2C. Average life of halogen bulbs : 4,000 Hrs. Power Requirement : 220V AC. Room Height : 9 ft



Samco Medical

Website: www.samcomedical.com,