



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ceiling Shadowless OT Light with Dia 504mm x 2 Domes,
Operated by Remote Control

Product Code :

SAMCOMEDICAL4007A029



Description :

Ceiling Shadowless OT Light with Dia 504mm x 2 Domes, Operated by Remote Control Key Structure A new superior quality O.T. Light executed by a Remote Control. This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated automatically in case the main light source brakes down. Sterilizable handle has been provided in the middle of the dome. Inbuilt is an Imported Dichoric Glass Reflector, which increases the Lux up to 1,80,000 (i.e. 10%). LED Display that helps the user to adjust the Intensity of Light from 0 to 100% & mode system. Technical specifications Dome head : 504 mm x 2 Halogen Bulbs : 24V x 150W x 4 nos. Light Intensity at 1 mt : 1,80,000 Lux Height Adjustment : 978mm Dichoric Glass Reflector : 390mm Color Temperature : 4200 K i.e. 1/2 300 Temp rise in field : 2i.e. 1/2 - 3i.e. 1/2 C Control Unit : SMPS



Samco Medical

Website: www.samcomedical.com,