



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Lensometer

Product Code :

SAMCOMEDICAL240001A013



Description :

Manual Lensometer: LED illumination and cord-free battery power. The cool-running, energy-efficient LED light source provides a clear green illumination target, and virtually maintenance-free operation. Battery operated and completely cordless. To verify the correct prescription in a pair of eyeglasses, to properly orient and mark uncut lenses, and to confirm the correct mounting of lenses in spectacle frames. Lensmeters can also verify the power of contact lenses Lensometer, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,