



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Suction Apparatus Trolley Model

**Product Code :**  
SAMCOMEDICAL390001A003



**Description :**

Suction Apparatus Trolley Model Ideal & suitable for Hospitals, Nursing Homes and Clinics  
• Cabinet made of mild steel sheet with epoxy powder coated  
• 1/2 HP heavy duty motor. Oil immersed vacuum pump  
• Two wide mouth glass bottles of approx. 1600 ml. each  
• Synthetic air tight rubber lids. Overflow cut-off device  
• Non-collapsible PVC Tubing  
• Vacuum-700 +/- 10mm Hg  
• Trolley mounted on friction free wheels  
Suction Apparatus Trolley Model, Suction Apparatus Exporters, Suction Apparatus Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),