



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vertical Sterilizer

Product Code :
SAMCOMEDICAL410001A010



Description :

Vertical Sterilizer We are among the reputed organizations, engaged in offering our clients with Vertical Sterilizer. All our products are widely demanded by several hospitals, nursing homes and clinics. These products are also checked by expert quality controllers on various parameters to ensure flawlessness. We also provide customization as per the specifications laid down by the customers. Features: •½ Corrosion & abrasion resistance •½ Excellent performance •½ Easy handling Specifications: Steam Working Pressure 1.26 kgf/cm² (2.2 kgf/cm² in case High speed sterilizer Steam Working Temperature 121oC (134o in case of High speed sterilizer) Hyd. Tested at Jacket-Twice the working pressure. Chamber-One & half time the working pressure. Operating Voltage 400/440, 3 phase, AC supply, 50c/s. Sterilization Period 20 to 25 minutes(in case of High speed sterilizer 5 to 7 minutes). Steam Exhaust 5 to 7 minutes(in case of High speed sterilizer 1 minute) Vertical Sterilizer, Infant Care Products Exporters, Infant Care Products Manufacturers



Samco Medical

Website: www.samcomedical.com,