



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Autoclaves Pressure Steam Sterilizer (Cooker Type)

Product Code :

SAMCOMEDICAL5006A001

Description :

Autoclaves Pressure Steam Sterilizer (Cooker Type)

Technical Specification :

Autoclaves Pressure Steam Sterilizer: Sizes Available : 300mm x 300 mm (12in x 12in) Made of Heavy Duty Aluminium, Double Safety Valve. Fitted with 2.5in Special Gauge. Steam Release Valve, Water Outlet Assembly, Working Pressure 15 lb Electric Load 2 KW Autoclaves Pressure Steam Sterilizer (Cooker Type), Autoclave & Sterilizers Exporters, Autoclave & Sterilizers Manufacturers



Samco Medical

Website: www.samcomedical.com,