



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Variable Double Beam Spectrophotometer Device

**Product Code :**  
SAMCOMEDICAL650001A002



**Description :**

Variable Double Beam Spectrophotometer Device Wavelength Range 190-1100 nm (1200 Lines/mm grating)  
Mode Basic/Quantitative/Wavelength Scan/DNA Protein Test/Kinetics Scanning Speed Fast / Medium / Slow  
Wavelength Accuracy +0.1 nm Bandwidth 0.5, 1,2,4,5 Wavelength Reproducibility 0.1 nm Photometric Display  
Range 0-200 % T, -0.3 - 3.0 A, 0-9999 C Photometric Accuracy + 0.3 % T Photometric Repeatability + 0.2 % T  
Stability 0.001A/h @ 500 nm Baseline Flatness + 0.001 A Noise +0.001A Stray Light < 0.05% T @ 220nm, 360  
nm) Data Output Port USB Printer Port Parallel Port Display Graphic LCD (320\*240 Dots) Lamps Deuterium  
Lamp & Tungsten Halogen Lamp Detector Silicon Photodiode Dimensions(L\*W\*H) 630 x 430 x 206 mm Net  
Weight 26 Kg Variable Double Beam Spectrophotometer Device, Double Beam UV-VIS Spectrophotometer  
exporters, Double Beam UV-VIS Spectrophotometer manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),