



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dental Compressor

Product Code :
SAMCOMEDICAL1400001A003



Description :

Dental Compressor 1/2 HP Dental compressor seal unit fitted with safety valve for air release excess sorted, auto cut off system 1/2 pressure gauge and drain valve air rotor air compressor mounted on high pressure tank entirely powder coated paint 1/2 Power 1/2 HP Amperes 4.5 Flow 40 to 50 Liter per minute 1/2 Speed 1000 rpm 1/2 Noise 50 to 60 dB (A) 1/2 Power supply 220V, 50 Hz 1/2 Tank 25L 1/2 Pressure 350 MMHG 1/2 Weight 35 to 37kgs



Samco Medical

Website: www.samcomedical.com,